



# Q-FLOW™ 1

## Compact examination light




- Single hand use
- Excellent color rendering capability
- Detachable sterilizable handle
- Easy to clean



# Q-FLOW™ 1

Q-Flow 1 is at it's best in any patient treatment area, such as: ward, emergency, examination, recovery and ICU.



Mobile with 2 antistatic castors



Ceiling mounted



Wall mounted

## Technical specifications

Max. illumination intensity at 1.0m* (lux)	50 000
Colour temperature (K)	4 400
Colour rendering index (Ra)	98
Red colour rendering index (R9)	93
“Skin colour” rendering index (R13)	98
Number of LEDs (pcs)	9
Power consumption (W)	17
Dimming (%)	25 – 100
Average life time of LED (h)	> 60 000
Working distance (mm)	500 – 1500
Depth of illumination 60 % (mm)	1 100
Light field diameter (mm)	240
Sterilisable handle	√
Size of light head (mm)	Length 400 Ø 260
Weight of light head (kg)	2,6
Degree of protection (lamp head)	IP44
Manufacturer	Merivaara Corp.

√ = included

\*) Tolerance ±10%



Easy to adjust and control with one hand.



**Hoyland Medical**

Service and quality to the medical profession

Unit 7, 22 Palmer Place

Murarie QLD 4172

T: 07 3390 7011

E: [admin@hoylandmedical.com.au](mailto:admin@hoylandmedical.com.au)