


Q-FLOW™ 2

Compact surgical light for
minor procedures



- High illumination power
- Excellent color rendering capability
- Detachable sterilizable handle
- Easy to clean

Q-FLOW™ 2

Q-Flow 2 is ideal medical light for minor surgical procedures in hospitals and private clinics.



Wall mounted



Ceiling mounted



Mobile with 2 antistatic castors

Technical specifications

| | |
|--|------------------|
| Max. illumination intensity at 1.0m* (lux) | 100 000 |
| Colour temperature (K) | 4 400 |
| Colour rendering index (Ra) | 98 |
| Red colour rendering index (R9) | 97 |
| “Skin colour” rendering index (R13) | 98 |
| Number of LEDs (pcs) | 18 |
| Power consumption (W) | 32 |
| Dimming (%) | 25 – 100 |
| Average life time of LED (h) | > 60 000 |
| Working distance (mm) | 500 – 1500 |
| Depth of illumination 60 % (mm) | 900 |
| Light field diameter (mm) | 208 – 220 |
| Sterilisable handle | √ |
| Size of light head (mm) | Length 400 Ø 260 |
| Weight of light head (kg) | 2,7 |
| Degree of protection (lamp head) | IP44 |
| Manufacturer | Merivaara Corp. |

√ = included

*) Tolerance ±10%



Single hand use together with double yoke option enables smooth positioning of the lamp head.



Unit 7, 22 Palmer Place

Murarrie QLD 4172

T: 07 3390 7011

E: admin@hoylandmedical.com.au