

Sonicaid FD1 / FD3 Dopplers

Rate display fetal Dopplers with fixed 2 or 3 MHz waterproof probe.



Large Monochrome Display

- BIG Numbers fetal heart rate display
- Heartbeat and battery indicators



Audio Output

- Headphone connection socket




High Sensitivity Probes

- Calibrated and fine-tuned to provide excellent performance over a wide gestational age & BMI range



Controls

- Combined Power ON/OFF/Volume control for ease of use




Probe Parking

- Protection when probe not in use



Hear the Difference

- High quality audio with analogue noise reduction



Integrated Belt/Pocket Clip

- For convenient carrying and storage

Accessories

The Sonicaid FD1/FD3 Dopplers are supplied with:

- A choice of 2 or 3MHz fixed waterproof probe
- 2 x 20g Sachets of ultrasound gel
- 2 x AA Alkaline batteries
- IFU

Optional accessories include:

- ACC34 - Soft carry pouch
- ACC52-2 - Mobile Doppler stand
- ACC3 - Ultrasonic gel (box 12x250ml bottles)
- ACC24 - Ultrasonic gel (box 12x60g tubes)
- ACC21 - Headset