




■ ■ **039 BASKET RACKS**




Wire diameter 6.0mm and 2.5mm

<u>New Item Code</u>	<u>Size in mm</u>	<u>Old Item Code</u>
CIL-039 001 001	600 x 400 x 15	CIL-006 500 015



Wire diameter 6.0mm and 3.0mm

<u>New Item Code</u>	<u>Size in mm</u>	<u>Old Item Code</u>
CIL- 039 002 001	600 x 400 x 44	CIL- 006 525 044



Wire diameter 6.0mm and 3.0mm

<u>New Item Code</u>	<u>Size in mm</u>	<u>Old Item Code</u>
CIL- 039 003 001	600 x 400 x 44	CIL- 006 550 044

- *The sterilization baskets are placed on racks during sterilization.*