



006 BASKETS SPRI

Double Frame with folding handle



- Upper and lower frame wire diameter 5.0 mm.
- Loops wire diameter 3.5 mm.



(With folding handle)

New Item Code	Size in mm	Old Item Code
CIL-006 001 001	585 x 395 x 100	CIL-006 001 100
CIL-006 001 002	585 x 395 x 195	CIL-006 001 195
CIL-006 002 001	600 x 300 x 97	CIL-006 003 097
CIL-006 002 002	600 x 300 x 138	CIL-006 003 138
CIL-006 002 003	600 x 300 x 150	
CIL-006 002 004	600 x 300 x 260	CIL-006 003 260
CIL-006 002 005	600 x 300 x 300	


New Item Code	Size in mm	Old Item Code
CIL-006 003 001	600 x 400 x 100	CIL-006 051 100

Double Frame with large folding handle



- Space for placing nomenclature plate.
- Upper and lower frame wire diameter 5.0 mm.
- Loops wire diameter 3.5 mm.

Double Frame




- Upper and lower frame wire diameter 5.0 mm.
- Loops wire diameter 3.5 mm.

New Item Code	Size in mm	Old Item Code
CIL-006 004 001	600 x 400 x 100	CIL-006 101 100
CIL-006 004 002	600 x 400 x 130	CIL-006 102 130
CIL-006 004 003	600 x 400 x 150	CIL-006 103 150
CIL-006 004 004	600 x 400 x 180	CIL-006 104 180

New Item Code	Size in mm	Old Item Code
CIL-006 005 001	600 x 400 x 100	CIL-007 001 100
CIL-006 005 002	600 x 400 x 130	CIL-007 001 130
CIL-006 005 003	600 x 400 x 150	CIL-007 001 150
CIL-006 005 004	600 x 400 x 180	CIL-007 001 180
CIL-006 005 005	600 x 400 x 200	CIL-007 001 200
CIL-006 006 001	600 x 300 x 50	CIL-007 002 50
CIL-006 006 002	600 x 300 x 100	CIL-007 002 100
CIL-006 006 003	600 x 300 x 130	CIL-007 002 130
CIL-006 006 004	600 x 300 x 150	CIL-007 002 150
CIL-006 006 005	600 x 300 x 180	CIL-007 002 180
CIL-006 006 006	600 x 300 x 200	CIL-007 002 200

Single Frame



- Nestable/Stackable with change of position of placing of baskets on each other
- Upper and lower frame wire diameter 6.0 mm
- Loops wire diameter 4.0 mm

Common Features

- The SPRI baskets are both stackable and nestable
- Empty baskets occupy minimum storage space while nesting position
- Stackable when placed on one another in opposite position
- Electrolyte polish and smooth finish