





**004 FINE MESH BASKETS WITH LOCKING CLIP**

<u>New Item Code</u>	<u>Size in mm</u>	<u>Old Item Code</u>
CIL-004 001 001	40 x 40 x 20	
CIL-004 001 002	40 x 40 x 30	
CIL-004 002 001	70 x 70 x 30	
CIL-004 003 001	80 x 40 x 20	CIL-004 001 020
CIL-004 003 002	80 x 40 x 30	CIL-004 002 030
CIL-004 004 001	80 x 80 x 34	CIL-004 003 034
CIL-004 004 002	80 x 80 x 40	CIL-004 004 040
CIL-004 004 003	80 x 80 x 50	
CIL-004 005 001	90 x 90 x 50	
CIL-004 006 001	100 x 100 x 50	
CIL-004 007 001	105 x 70 x 25	CIL-004 005 025
CIL-004 008 001	150 x 50 x 35	CIL-004 006 035
CIL-004 009 001	190 x 85 x 100	CIL-004 007 100
CIL-004 010 001	200 x 100 x 35	CIL-004 008 035
CIL-004 011 001	220 x 140 x 50	CIL-004 009 050
CIL-004 012 001	235 x 155 x 40	CIL-004 010 040
CIL-004 013 001	240 x 90 x 40	CIL-004 011 040
CIL-004 014 001	240 x 120 x 35	CIL-004 012 035
CIL-004 015 001	275 x 175 x 35	CIL-004 013 035
CIL-004 016 001	275 x 200 x 35	CIL-004 014 035



*Stainless steel fine mesh baskets are used for micro medical and surgical instruments.*

- Top and lid frame wire diameter 3.5 mm and 4.0 mm (according to the size of the baskets).
- Fine mesh with opening square 1.6 mm and wire diameter 0.7 mm.
- Wire locking clip.
- Electrolyte polish.
- Custom sizes are offered.



**034 ULTRA FINE BASKET**

*Ultra fine baskets are used for delicate micro surgical and laboratory instruments for sterilization.*

- Top and lid frame wire diameter 3.0 mm
- Ultra fine mesh with opening square 0.5 mm and wire diameter 0.3 mm.
- Wire locking clip.
- Electrolyte polish.



<u>New Item Code</u>	<u>Size in mm</u>	<u>Old Item Code</u>
CIL-034 001 001	80 x 40 x 30	CIL-004 202 030
CIL-034 002 001	80 x 80 x 34	CIL-004 201 034
CIL-034 003 001	180 x 130 x 20	