




■ ■ **002 LID FLAT BASE BASKETS**



<u>New Item Code</u>	<u>Size in mm</u>	<u>Old Item Code</u>
CIL-002 001 501	120 x 250	CIL-002 001 020
CIL-002 002 501	240 x 250	CIL-002 002 020
CIL-002 003 501	255 x 245	CIL-002 003 020
CIL-002 004 501	405 x 250	CIL-002 004 020
CIL-002 005 501	405 x 255	CIL-002 005 020
CIL-002 006 501	480 x 250	CIL-002 006 020
CIL-002 007 501	485 x 255	CIL-002 007 020
CIL-002 008 501	480 x 360	CIL-002 008 020
CIL-002 009 501	480 x 480	CIL-002 009 020
CIL-002 010 501	515 x 245	CIL-002 010 020
CIL-002 011 501	540 x 250	CIL-002 011 020



The flat base baskets are used as sterilizer accessories for medical & surgical products in laboratory sterilizers

- *Separate lids are available.*
- *Available in various dimensions.*
- *Burr free side walls with perforated sheet*
- *Base in **wire net** with opening square 6.0 mm and wire diameter 1.2 mm*
- *Base also available in **wire mesh** with opening square 4.2 mm and wire diameter 0.9 mm*
- *Inner mounted drop handles*
- *Stay Handle*
- *Smooth finish electrolyte polish*
- *Custom sizes are offered.*