

Dopplex[®] SD2 Doppler

Provides bi-directional blood flow visualisation




Large Monochrome Display

- Large monochrome display screen with blood flow direction indicator to easily identify forward and reverse flow
- Clear probe frequency identification



Audio Output

- Flow separated stereo headphone output



High Sensitivity Interchangeable Probes

- All probes have been calibrated and fine-tuned for optimum use




Controls

- Combined Power ON/OFF/Volume control for ease of use




Hear the Difference

- Clear audio sounds of blood flow



Probe Parking

- Protection when probe not in use



Integrated Belt/Pocket Clip

- For convenient carrying and storage



Accessories

The Dopplex SD2 is supplied with:

- 2 x 20g Sachets of ultrasound gel
- 2 x AA Alkaline batteries
- IFU

Optional accessories include:

- ACC21 - Stereo headset
- ACC34 - Soft carry pouch
- ACC52-2 - Doppler stand
- ACC3 - Ultrasonic gel (box 12x250ml bottles)
- ACC24 - Ultrasonic gel (box 12x60g tubes)
- ACC188 - Blood pressure cuff & sphyg

Probes: (Sold separately)

- VP4 - 4MHz
- VP5 - 5MHz
- VP8 - 8MHz
- EZ8 - 8MHz Widebeam
- VP10 - 10MHz
- Surgical probes (refer to page 10)