




■ ■ **005 STERILE GOODS BASKETS**

Standard



<u>New Item Code</u>	<u>Size in mm</u>	<u>Quality</u>	<u>Old Item Code</u>
CIL-005 001 001	575 x 280 x 135	Standard	CIL-005 001 135
CIL-005 001 002	575 x 280 x 265	Standard	CIL-005 001 265
CIL-005 001 003	575 x 280 x 265	Front access	CIL-005 002 265
CIL-005 002 001	600 x 300 x 150	Standard	
CIL-005 002 002	600 x 300 x 260	Front access	
CIL-005 002 003	600 x 300 x 300	Standard	
CIL-005 002 004	600 x 300 x 300	Front access	
CIL-005 003 001	600 x 400 x 150	Standard	
CIL-005 003 002	600 x 400 x 300	Standard	
CIL-005 003 003	600 x 400 x 300	Front access	

Front access



The standard baskets are used for all types of sterilizers. STE sterile goods baskets mean the baskets made with continuous wire which provide maximum protection to user from being injured. The baskets are stackable and can also be mounted on the cart or wall with hooks.

- Upper frame wire diameter 6.0 mm.
- Frame wire diameter 3.0 mm.
- Loops wire diameter 2.5 mm.
- Bottom frame wire diameter 5.0 mm.
- Electrolyte polish and burr free smooth finish.