





## 009 WIRE NET BASKET



<u>New Item Code</u>	<u>Size in mm</u>	<u>Old Item Code</u>
CIL-009 001 001	250 x 170 x 70	CIL-009 001 070
CIL-009 002 001	340 x 250 x 70	CIL-009 002 070
CIL-009 003 001	450 x 340 x 70	CIL-009 003 070

Wire net baskets are used for sterilization of instruments in a sterilizer.

- Bottom frame fabricated to allow space between instruments and working surface
- Separate top covers are available
- Stackable
- Top and bottom frame wire diameter 3.5 mm
- Grid made from wire diameter 1.6 mm having opening square 6.4 mm
- Burr free smooth finish
- Electrolyte polish



## 009 LID WIRE NET

<u>New Item Code</u>	<u>Size in mm</u>	<u>Old Item Code</u>
CIL-009 001 501	250 x 170	CIL-009 004 04
CIL-009 002 501	340 x 250	CIL-009 005 04
CIL-009 003 501	450 x 340	CIL-009 006 04