



Spectra Range

The Spectra range is a cart system that combines mobility and drawer flexibility into a lightweight colour coded package that has applications throughout all healthcare environments. The unique drawer system is based on removable 60 x 40 tray systems that offer different depths and the option of integrated drawer dividers. This allows flexibility with the carts applications, particularly in the ability to restock the carts with transfer tray systems.

The carts shelves are one piece moulded polymer which enables easy cleaning and eliminates infection control issues due to the rounded internal mould design, whilst integrating a raised lip to three sides to prevent items falling off.

The Spectra range offers unlimited drawer configurations and a range of colourful drawer fronts.

Features:

- Anti bacterial polymer shelf
- Spectra colour drawer options
- Module 60 x 40 trays or baskets
- 125mm Castors; 2 Brake

Optional:

- Interlocking tray dividers
- Ergonomic push handle




AXS0010

1 Drawer Trolley
 700 x 500 x 1000H mm


- 1 x 100mm deep tray
- Work surface area 620 x 460mm



Tray 60 x 40 x 5 cm



Tray 60 x 40 x 10 cm



Tray 60 x 40 x 20 cm



See page 2 for 60 x 40 Trays

Spectra Colour Options



Red

Yellow

Green

Blue



AXS0020
2 Drawer Trolley
700 x 500 x 1000H mm

- 2 x 100mm deep trays
- Work surface area 620 x 460mm

AXS0030
Multi Drawer Trolley
700 x 500 x 1000H mm

- 1 x 50mm deep tray
- 2 x 100mm deep trays
- 2 x 200mm deep trays
- Work surface area 620 x 460mm

