

BD 48

Electronic Blood Drawing Chair, 4 Motors



BD 48 Blood Transfusion Chair designed for blood transfusion, hemodialysis, eye operations, skin care, body care, epilation, ozone therapy, oncology and other medical esthetic operations.

- 4 motors
- Electronic adjustments of backrest, footrest, height, Trendelenburg and reverse Trendelenburg functions
- Hand control unit
- Auto-contour (simultaneous adjustments)
- Tilt proof
- Water-proof, easy cleanable synthetic leather upholstered seat and back surface
- ABS covered sub-frame and outer sides of chair
- Synthetic leather upholstered armrests
- Lockable castors
- Adjustable Foot end
- Battery backup (optional)

TECHNICAL DETAILS	
Width with armrests:	97 cm
Length:	181 cm
Height:	64-84 cm
Seating Area Width:	56 cm
Backrest Angle:	0°-75°
Footrest Angle:	0°-88°
Trendelenburg Angle:	0°-13°
Reverse Trendelenburg Angle:	0°-5°
Safe Working Capacity:	210 kg
Weight:	100kg