



Polyurethane wrap to preserve the body's own warmth and moisture



Premature and full-term babies can lose a disproportionate amount of heat and moisture due to their immature skin.

neoTHERM® is a very simple but effective way of preventing heat and moisture loss by wrapping the baby in a wafer-thin protective layer. The integrated fleece gently absorbs the moisture and thus offers additional skin protection.

neoTHERM® can be used in a number of ways: as a standard in primary care or part of procedures such as delayed cord clamping. During intrahospital transport and examinations such as MRI or for simple infection protection.

A key benefit of the neoTHERM® is that the mother can continue to have intensive skin-to-skin contact with her little treasure whilst using it – 'kangaroo' care is possible without restriction.

#### Product specifications

- Single patient use
- Available in three sizes
- Surface material: polyurethane film
- Double insert: removable fleece
- Medical device class 1

#### Innovations of the neoTHERM®

- Soft and supple
- Skin-to-skin contact with mother possible
- Individual adaptation to patient size
- Easy handling
- Unrestricted patient observation
- Can be used in the incubator

#### Ordering information

PRODUCT DESCRIPTION	AMOUNT	ORD.-NO.:
neoTHERM® wrap for premature and full-term babies for first aid (anatomical design with hood), size S (< 1.000 g)	10 pcs.	TC9280210
neoTHERM® wrap for premature and full-term babies for first aid (anatomical design with hood), size M (1.000 g - 2.000 g)	10 pcs.	TC9280110
neoTHERM® wrap for premature and full-term babies for first aid (anatomical design with hood), size L (> 2.000 g)	10 pcs.	TC9280010